


# *Student Experiments*

Manual

# VIBRATIONS AND WAVES

P9160-4S



# INDEX

## 1. VIBRATIONS

- SWS 1.1.1 Oscillation period of a simple pendulum
- SWS 1.1.2 Oscillation period of a coil spring pendulum
- SWS 1.1.3 Oscillation period of a flat spring
- SWS 1.2 Path time chart of harmonic oscillation
- SWS 1.3 Measuring acceleration due to gravity
- SWS 1.4.1 Resonance of a simple pendulum
- SWS 1.4.2 Resonance of a coil spring pendulum
- SWS 1.4.3 Resonance of a flat spring
- SWS 1.5 Principle of a resonant vibrating-reed frequency meter
- SWS 1.6 Dynamic measuring of a spring constant

## 2. WAVES


- SWS 2.1 Stationary transversal wave
- SWS 2.2 Stationary longitudinal wave
- SWS 2.3 Reflection of waves at their fixed and loose ends

# PATH TIME CHART OF HARMONIC OSCILLATION

SWS 1.2

## Required Kit:


P9901-4A Rail stand material  
P9901-4B Mechanics 1  
P9901-4S Vibrations and Waves



## Material:

1x Stand rail, 300mm  
1x Table clamp  
1x Sliding saddle for screens, springs and pointers  
1x Holder for pencil  
2x Slotted weights, 50g

1x Measuring tape  
1x Threaded rod with butterfly nut  
1x Flat spring steel  
Pencil, Sheet of paper A4



# PATH TIME CHART OF HARMONIC OSZILLATION

SWS 1.2

This temporal course of oscillation can be seen by introducing a time axis. In this experiment the time axis is produced by uniform movement of a strip of paper during the oscillation of the flat spring.


**Preparation:** Arrange according to the illustration. The table clamp with stand rail is fixed to the edge of the table. The sliding saddle for scales is placed on the stand rail. The flat spring is clamped in this sliding saddle. At the end of the flat spring the holder for the pencil is mounted in such a way that the drill hole of the holder is exactly above the drill hole of the flat spring. The threaded rod is pushed through this drill hole. One 50 g slotted mass is put on the left and one on the right hand side of the flat spring and fixed by means of the butterfly nut. The pencil is inserted in the holder. The flat spring is clamped in such a way that the distance of the pencil from the edge of the sliding saddle with setscrew is 30 cm.

A sheet of paper is placed beneath the flat spring. The pencil is to touch the sheet in its center.

**1. Experiment:** The spring is deviated about 4 cm. It starts oscillating. During the oscillation the paper is pulled as uniformly as possible (normal to the direction of oscillation) underneath the pencil.

**2. Experiment:** The spring is clamped in such a way that the distance of the pencil to the sliding saddle is half as long as in the first experiment (15 cm). The sheet of paper is placed beneath the spring in such a way that the recording of the oscillation can start at the same place as the recording of the first oscillation. The flat spring is set in oscillation. The paper is moved along underneath the pencil at about the same speed as in the first experiment.

**Conclusion:** The temporal course of the oscillation of a flat spring is sinusoidal oscillation. The frequency of the oscillation depends on the length of the oscillating flat spring. It is the higher the shorter the length of the flat spring is.



# *Student Experiments*

© Fruhmann GmbH  
NTL Manufacturer & Wholesaler

Werner von Siemensstraße 1  
A - 7343 Neutal  
Austria

[www.ntl.at](http://www.ntl.at)