


Student Experiments

Manual

ALTERNATIVE ENERGY- CONVERSION

P9160-4W



INDEX

1. Combustion engines

- AES 1.1 Combustion process in a gasoline engine
- AES 1.2 Combustion process in a diesel engine

2. Energy made perceptible

- AES 2.1 Energy demand of a light bulb
- AES 2.1.1 Energy demand of a light bulb - quantitatively
- AES 2.1.2 Energy demand of home appliances
- AES 2.2 Energy demand of a small electric car
- AES 2.2.1 Additional experiments with the electric car
- AES 2.2.2 Electric car - practical cross reference
- AES 2.3 Energy demand of a propeller
- AES 2.3.1 Additional experiments with the propeller

3. Wind power

- AES 3.1 Wind power station
- AES 3.1.1 Wind power - continuative experiments

4. Water power

- AES 4.1 The hdyropower station
- AES 4.1.1 Water power - continuative experiments

5. Thermal power (caloric force)

- AES 5.1 Principle of thermal power plants
- AES 5.1.1 Thermal power - continuative experiments

6. Sun power - Photovoltaics

- AES 6.1 Voltage metering and incidence angle - qualitatively
- AES 6.2 Series connection of solar cells
- AES 6.3 Parallel connection of solar cells
- AES 6.4 Motor with propeller driven by a solar cell
- AES 6.5 Electric car driven by a solar cell
- AES 6.6 Angle of incidence and voltage

7. Energy storage

- AES 7.1 Quickly charging the short-time energy storage
- AES 7.2 The short-time energy storage as energy source
- AES 7.3 Charging facilities of the short-time energy storage



PRINCIPLE OF THERMAL POWER PLANTS

AES 5.1


You might know coal-fired power stations, gas- or oil-fired power plants. Such power plants are known as kaloric power plants or thermal power plants. But where does the electrical energy comes from in these kaloric power plants?

Experiment:


- Place the turbine in housing on the engine/generator brick
- Insert the metal shaft of the engine/generator in the turbine axis
- The 4 plastic bolts of the turbine housing have to stick in the opening of the engine/generator-brick. Stick both bricks together.
- Press the turbine with the black push button on the generator axis
- Check the smooth running of the turbine with your finger
- Screw the lamp 1.5V into the socket of the brick
- Connect both bricks with cables (red and blue, 75cm)



- Place the air pump on the ground next to the table
- Connect the flexible hose with the "pressure"-side of air pump
- Adjust the nozzle of the pressure hose to the turbine housing and onto the center of the "blade surfaces"



- Hold the turbine housing
- An other student assists by pumping



Result:

The air pressure drives the turbine, the lamp starts to glow.

PRINCIPLE OF THERMAL POWER PLANTS

AES 5.1

Physical explanation of the energy conversion:

Mechanical energy (rotation of the turbine) is converted into electrical energy (generator). This electrical energy is then converted into light- and heat-energy (lamp).

Cross-reference to practise:

In a thermal power plant water gets heated-up in a container through burning fuel. This causes the water to evaporate, high pressure occurs in the container. The hot steam escapes from this high pressure through a nozzle and drives the turbine.


Try to describe a thermal power station (operated by coal, gas or oil) schematically. Draw the vapor pressure nozzle and a turbine too.

Advice:

As the operation of the turbine through hot gas steam is difficult to manage in a laboratory experiment (also this would be very dangerous because of the high temperatures, we illustrate the principle by using an air pump.

Additional explanation of nuclear power plants:

Even in nuclear power plants a liquid gets heated-up. This heat-energy gets then converted by a turbine into electrical energy. The difference to "conventional" thermal power plants can be found in the combustion process.



Student Experiments

© Fruhmann GmbH
NTL Manufacturer & Wholesaler

Werner von Siemensstraße 1
A - 7343 Neutal
Austria

www.ntl.at