

VSC[®] MSFI-2

Magnetic Security Feature Imager

Compact USB device for the examination and advanced verification of magnetic security features on banknotes and other secure documents.

MSFI-2 is a high-quality device for the visualisation of the hard and soft magnetic inks that are commonly used in the printing of secure documents including banknotes and cheques.

Manufactured from high quality components the MSFI-2 includes a detachable magnetic excitation system, tuneable cross polariser for contrast adjustment, integral illumination source, and a single USB port for connection to a VSC workstation.



MSFI-2

Order # V8/MSFI/2

Comprises:

- MSFI 2 Main Unit
- MSFI 2 Docking Tray
- Magnetic Excitation System
- USB 2.0 cable
- System Software